



TB

BIPOLAR TEMPORARY
PACING LEADS

LATEX-FREE



TB Bipolar Temporary Pacing Lead:

TB is a bipolar temporary pacing lead used for temporary ventricular or atrial pacing in conjunction with a temporary pacemaker. Oscor's TB Series are available in a variety of curve configurations, and french sizes.

TB Features:

- Lead body diameter of 4 F, 5 F, or 6 F
- Usable lead length 110 cm
- Polymer radiopaque lead material
- 10 mm electrode spacing
- Stainless steel electrode material
- Electrical resistance of < 100 Ω
- Atraumatic tip configuration
- Fixed connectors available with 2 mm pins, shrouded or unshrouded

Model # Kit	French Size	Pin Configuration	ATAR™ Cable	Suggested Introducer
<i>(One per box)</i>				
TB—Straight Curve (Ventricular Pacing)				
TB 4 F S	● 4 F	Shrouded	ATAR™ T	4 F
TB 5 F S	● 5 F	Shrouded	ATAR™ T	5 F
TB 6 F S	● 6 F	Shrouded	ATAR™ T	6 F
TB 4 F U	● 4 F	Unshrouded	ATAR™ V	4 F
TB 5 F U	● 5 F	Unshrouded	ATAR™ V	5 F
TB 6 F U	● 6 F	Unshrouded	ATAR™ V	6 F
TB V—60° Curve (Ventricular Pacing)				
TBV 4 F S	● 4 F	Shrouded	ATAR™ T	4 F
TBV 5 F S	● 5 F	Shrouded	ATAR™ T	5 F
TBV 6 F S	● 6 F	Shrouded	ATAR™ T	6 F
TBV 4 F U	● 4 F	Unshrouded	ATAR™ V	4 F
TBV 5 F U	● 5 F	Unshrouded	ATAR™ V	5 F
TBV 6 F U	● 6 F	Unshrouded	ATAR™ V	6 F
TB J—Atrial-J (Atrial Pacing)				
TBJ 4 F S	● 4 F	Shrouded	ATAR™ T	4 F
TBJ 5 F S	● 5 F	Shrouded	ATAR™ T	5 F
TBJ 6 F S	● 6 F	Shrouded	ATAR™ T	6 F
TBJ 4 F U	● 4 F	Unshrouded	ATAR™ A	4 F
TBJ 5 F U	● 5 F	Unshrouded	ATAR™ A	5 F
TBJ 6 F U	● 6 F	Unshrouded	ATAR™ A	6 F
TB RH—Right Heart Curve (Ventricular Pacing)				
TBRH 4 F S	● 4 F	Shrouded	ATAR™ T	4 F
TBRH 5 F S	● 5 F	Shrouded	ATAR™ T	5 F
TBRH 6 F S	● 6 F	Shrouded	ATAR™ T	6 F
TBRH 4 F U	● 4 F	Unshrouded	ATAR™ A	4 F
TBRH 5 F U	● 5 F	Unshrouded	ATAR™ A	5 F
TBRH 6 F U	● 6 F	Unshrouded	ATAR™ A	6 F

Leads supplied sterile. For single use only.

Standard lengths are noted above. Other lengths available upon request

For additional information e-mail us at sales@oscor.com.



P: 727.937.2511
F: 727.934.9835
E: sales@oscor.com
www.oscor.com

CAUTION: Federal Law restricts this device to sale by or on the order of a physician (U.S.A). Please refer to the Instructions for Use for a complete listing of the indications, contraindications, precautions and warnings, where applicable.

